BIOLOGIA MOLECOLARE

NLM 7119

VER 03/06/2014

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REF AC004/24

24 TEST

CND W01050203

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BIBLIOGRAFIA

1. EASL Clinical Practice Guidelines: Management of hepatitis C virus infection. *Journal of Hepatology* 2011 vol.55: 245-264.
2. Nizar N. Zein. Clinical Significance of Hepatitis C Virus Genotypes. *Clinical Microbiology Reviews*, Apr. 2000, p.223-235.
3. Verbeeck J., et al. Use of a Commercially Available Line Probe Assay for Genotyping of Hepatitis C Virus 5a Strains. *Journal of Clinical Microbiology*, Dec. 2005, p.6117-6119
4. Simmonds et al. Consensus proposals for a unified system of nomenclature of hepatitis C virus genotypes. *Hepatology*, vol. 42, n° 4, 2005.
5. Hino K, Yamaguchi Y, Fujiwara D, Katoh Y, Korenaga M, Okazak, Ohuda M, Okita K. Hepatitis C virus quasispecies and response to interferon therapy in patients with chronic hepatitis C: a prospective study. *J Viral Hepat* 2000; 7:36-42.
6. Khan MH, Farrell GC, Byth K, Lin R, Weltman M, George J et al. Which patients with hepatitis C develop liver complications. *Hepatology* 2000; 31:513-520.
7. Guideline on clinical evaluation of medicinal products for the treatment of chronic hepatitis C. *European Medicines Agency*. Jan. 2011
8. Hadziyannis SJ. Why and how treat chronic hepatitis C. *Can J Gastroenterol* 2000; 14 suppl B:45B-48B
9. Griso D., Rivanera D, Lilli D., Macri A., Biolcati G., Piunno M., Misogano N., Lozzi M.A., Mancini C.. Role of viral hepatitis infection in Italian patients with sporadic and familial porphyria cutanea tarda. Millenium meeting on porphyrins and porphyrias 2000. Institut Pasteur Paris – France. September 10-13, 2000.
10. P.T. Hraber et al. Comparative analysis of hepatitis C virus phylogenies from coding and non-coding regions: the 5' untranslated region (UTR) fails to classify subtypes. *Virology Journal* 2006, 3:103
11. M. A. Ansari et al. HCV-Core Region: Its Significance in HCV-Genotyping and Type Dependent Genomic Expression. *Macedonian Journal of Medical Sciences*. 2012 Mar 15; 5(1):30-39.

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24 TEST

GMDN 30742

AMPLIFICATION: AA910/48 or AA896/48
EXTRACTION: AA1318 NOT INCLUDED

GEN-C 2.0

Reverse Hybridization Strip Assay

CE

0459

IVD

NUCLEAR LASER MEDICINE S.r.l.
HEAD OFFICE: Viale delle Industrie, 3 – 20090 SETTALA MI (Italy)
Tel (+39) 02/95. 24. 51 - Fax (+39) 02/95. 24. 52. 37
WEB: www.nlm.it – E-MAIL: info@nlm.it

⇒ INTENDED USE

Gen-C 2.0 is an *in vitro* line probe assay for the genotyping of hepatitis C virus. The assay discriminates between HCV genotypes on the basis of variations in the 5'-UTR and Core regions. The assay works with amplicons obtained with NLM code AA910/48 or AA896/48. For optimal results it is strictly recommended to use samples with viral load between 5×10^3 and 8×10^6 (if necessary dilute high viral load plasma before the RNA extraction).

On the basis of the *in silico* sequence analysis, the device is able to distinguish genotypes 1, 2, 3, 4 and 6, and subtypes 1a, 1b, 2a/2c, 2b, 3a, 3b, 3c, 3k, 4a, 4b, 4c/4d, 4e, 4f, 4h, 5a, 6a/6b, 6g,f,q, 6m and 7a.

The result of this assay should be intended as an aid in selecting type and length of antiviral therapy for patients with chronic HCV infection. Therefore the device is not intended for use as a screening test for the presence of HCV RNA or as a diagnostic test to confirm the presence of HCV infection.

⇒ INTRODUCTION

A very high number of HCV-infected patients develop chronic hepatitis which often results in liver cirrhosis and hepatocellular carcinoma (1). HCV genome is a linear single strand RNA consisting of about 9400 nucleotides. The Core, Envelope and non-structural domains form the coding region of the genome, flanked in 5' and 3' by highly conserved untranslated regions (UTR).

The genetic heterogeneity of coding and untranslated regions determines the classification of HCV into different genotypes (2,3). Current nomenclature allows the classification in: genotypes, which have nucleotides difference of about 31-33%, and subtypes, that exhibit nucleotide differences of about 20-25%. This nomenclature system is widely used by the international scientific community (4) and define the following genotypes: 1a to 1c, 2a to 2d, 3a to 3k (previously defined as 10a), 4a to 4k, 5a, 6a-q 7a. The 5'UTR region is characterized by variable regions among the different HCV genotypes that are useful for distinguishing types 1-5, 6a-b and 7a, but is not able to differentiate accurately subtypes 6g,f,q, 6m and subtypes 1a and 1b. Core region elements were added to better identify these types and to improve genotype 1 subtyping accuracy (10-11).

IFN- α treatment

Results from many clinical studies regarding interferon (IFN) treatment of chronic HCV indicated that different genotypes respond in a dissimilar way if treated with IFN- α . In particular genotype 1b HCV infection could not be efficiently treated with IFN- α (less than 10% of treated patients with a long-term response), while genotypes 1a, 2a, 2b, and 3a responded favorably (between 50 and 80% long-term response). It has been demonstrated that genotype 1b infections proceed much faster to severe forms of chronic hepatitis, cirrhosis, and hepatocellular carcinoma (5,6). Treatment response is dependent by age, sex, and duration of disease (6). New drugs like Ribavirin together with IFN- α improved the number of therapeutic success (7-9). Discrimination of HCV subtypes is of great diagnostic importance since it helps to find more correct, more effective and personalized therapeutic protocols.

REFERENCES

1. EASL Clinical Practice Guidelines: Management of hepatitis C virus infection. *Journal of Hepatology* 2011 vol.55: 245-264.
2. Nizar N. Zein. Clinical Significance of Hepatitis C Virus Genotypes. *Clinical Microbiology Reviews*, Apr. 2000, p.223-235.
3. Verbeeck J., et al. Use of a Commercially Available Line Probe Assay for Genotyping of Hepatitis C Virus 5a Strains. *Journal of Clinical Microbiology*, Dec. 2005, p.6117-6119
4. Simmonds et al. Consensus proposals for a unified system of nomenclature of hepatitis C virus genotypes. *Hepatology*, vol. 42, n° 4, 2005.
5. Hino K, Yamaguchi Y, Fujiwara D, Katoh Y, Korenaga M, Okazaki, Ohuda M, Okita K. Hepatitis C virus quasispecies and response to interferon therapy in patients with chronic hepatitis C: a prospective study. *J Viral Hepat* 2000; 7:36-42.
6. Khan MH, Farrell GC, Byth K, Lin R, Weltman M, George J et al. Which patients with hepatitis C develop liver complications. *Hepatology* 2000; 31:513-520.
7. Guideline on clinical evaluation of medicinal products for the treatment of chronic hepatitis C. *European Medicines Agency*. Jan. 2011
8. Hadziyannis SJ. Why and how treat chronic hepatitis C. *Can J Gastroenterol* 2000; 14 suppl B:45B-48B
9. Griso D., Rivanera D, Lilli D., Macri A., Biolcati G., Piunno M., Misogano N., Lozzi M.A., Mancini C.. Role of viral hepatitis infection in Italian patients with sporadic and familial porphyria cutanea tarda. Millenium meeting on porphyrins and porphyrias 2000. Institut Pasteur Paris – France. September 10-13, 2000
10. P.T. Hrabec et al. Comparative analysis of hepatitis C virus phylogenies from coding and non-coding regions: the 5' untranslated region (UTR) fails to classify subtypes. *Virology Journal* 2006, 3:103
11. M. A. Ansari et al. HCV-Core Region: Its Significance in HCV-Genotyping and Type Dependent Genomic Expression. *Macedonian Journal of Medical Sciences*. 2012 Mar 15; 5(1):30-39.

⇒ PRINCIPLE OF THE PROCEDURE

The test is based on reverse-hybridization principle. Biotinylated amplicons, generated by RT-PCR of the 5'-UTR and Core regions of HCV RNA, are hybridized to specific probes that are bound to nitrocellulose strip by a poly-T tail.

Biotinylated hybrids are then detected using streptavidin bound to alkaline phosphatase; amplicons that are not complementary are washed out. Then the substrate reacts with the streptavidin-alkaline phosphatase complex forming a purple precipitate and colouring banding pattern on the strip.

PRODUCT COMPOSITION

REAGENTS (2-8°C)	NLM Code	Quantity	N° VIAL
Strips Nitrocellulose membranes coated with oligonucleotides	KC029	24	1
DNAT R36/38* Denaturation solution containing NaOH	KA107	1 ml	1
Hybridization Saline Solution containing preservatives	KA102	60 ml	1
Stringent Wash Solution Saline solution containing detergent and preservatives	KA573	175 ml	1
Wash Solution B Saline Solution containing preservatives	KA105	175ml	1
Conjugate Solution containing streptavidin-labeled with alkaline phosphatase, stabilizers and preservatives	KA574	60 ml	1
Color developer Solution containing 5-Bromo-4-chloro-3-indolyl-phosphate and 4-Nitroblue tetrazolium	AA564	60 ml	1
Plastic Trays	DC000	3	
Interpretation chart Sheet that provides banding patterns for identification of genotypes		1	
Transparent interpretation card		1	

*Irritating to eyes and skin

⇒ MODULARITY

Reagents are provided for max. 3 distinct runs. If modularity isn't respected reagents could not be enough for the number of tests declared.

STORAGE AND STABILITY

- All opened or unopened reagents are stable up to the expiry date indicated on the label when stored at the correct temperature.
- At the end of each assay store reagents at the correct temperature.
- The kit should be kept isolated from any source of contaminating DNA especially amplified products.
- Developed dry strips should be stored in the dark at room temperature (+15/30°C).
- The vial containing the DNAT Solution should be closed immediately after use.

PRECAUTIONS

- Only professional and opportunely trained personnel should use this kit. Handle this product according to established good laboratory practices and universal precautions; wear personal protective apparel.
- All disposable items (tips and tubes) used during amplification step have to be sterile. When handling DNA samples or amplification products, use aerosol-resistant pipette tips to avoid pipettes contamination.
- Discard used materials as bio hazardous waste.
- Plastic trays can be reused maximum twice, if washed with bleach and then with water at the end of detection.
- Do not eat, drink, smoke or apply cosmetics in areas where reagents or specimens are handled.
- If skin or mucous membrane exposure occurs, immediately wash the area with copious amount of water. Seek medical advice.
- Do not use components beyond the expiration date.
- Test strips stored at +2/+8°C in the dark are stable until the expiry date.
- Do not mix reagents from different lots.
- Reagents have to be preserved separated from possible contaminants (as DNA samples and amplification products).
- Pay attention to device modularity: if modularity isn't respected, reagents will not be enough to process all the tests.
- It is recommended to work in three separated areas:
 - Area 1: pre-PCR (samples handling and extraction).
 - Area 2: Master Mix preparation.
 - Area 3: post-PCR (amplification and detection).
- Risk phrases linked to some components:
 - **R36/38**: irritating to eyes and skin
- Don't use the device if the box is damaged; contact the supplier.
- **It is advisable to have constant and uniform laboratory temperature, avoid to place the instruments near heating/cooling sources that may compromise the correct working.**
- **Use the transparent interpretation card of the specific lot.**

AMPLICONS STORAGE

Store amplicons obtained with AA910/48 ver.2.0 or AA896/48 ver.2.0 at +2/+8°C or at -15/-25°C for later use.

MATERIALS REQUIRED BUT NOT PROVIDED

- Water bath with inclined lid, shaking platform (50-100 rpm) and adjustable temperature (46°C ± 0,5 °C) or **Tecan ProfiBlot 48 with software ver. 1.xx**
- Thermometer
- Deionised or distilled water
- Dedicated adjustable volume micropipettes set and tips
- Clean tweezers to handle strips
- Vacuum aspiration apparatus

PROCEDURE

⇒ MANUAL DETECTION

- Adjust the water level of the shaking bath to approx. 2/3 of the height of the Plastic Tray.
- Heat the water bath to 46°C and check the temperature (46 ± 0,5°C) using a calibrated thermometer.
- Pre-warm the Hybridization and Stringent Wash Solution to 46°C; all crystals should be completely dissolved.
- Allow DNAT, Conjugate, Wash Solution B and Color developer to reach room temperature.
- Take one strip for each sample using clean tweezers (touch strips with gloves only). Label strips with a pencil (do not use ballpoint pen, markers, etc).
Warning: do not allow the strips to dry during the whole procedure.

Hybridization (46°C)

- Pipette **10 µl** of **DNAT** into each lane of the Plastic Tray
Warning: do not use DNAT if is not blue.
- For each sample add **10 µl** of **amplification** product to the corresponding drop of DNAT. Mix thoroughly with a pipette (solution will remain blue).
- Incubate **5 min** at room temperature.
- Add **1 ml** of **Hybridization** (pre-warmed to 46°C) to each sample and mix gently: blue color disappears.
- Place strips into the respective lanes with marker lines facing up. Submerge completely.
- Incubate for **30 min** at **46°C** in the shaking water bath (about 50 rpm). Keep the water bath closed to avoid temperature variations.

Stringent Wash solution (46°C)

- After hybridization step remove the tray from the water bath; angle the tray upwards and aspirate the solution from each lane using a pipette or a vacuum aspiration apparatus, without touching the strip surface to avoid its damage.
- Add **1 ml** of **Stringent Wash solution** (pre-warmed at 46°C) to each lane and rinse briefly.
- Remove the solution and add **1 ml** of **Stringent Wash solution** (pre-warmed at 46°C).
- Incubate **15 min** at **46°C** in the shaking water bath (about 50 rpm).
Warning: put Stringent Wash solution back in the water bath during incubation
- Remove the solution and add **1 ml** of **Stringent Wash solution** (pre-warmed at 46°C).
- Incubate **15 min** at **46°C** in the shaking water bath (about 50 rpm).
- Remove the solution from each lane.

Color development (room temperature)

All the next steps have to be performed at room temperature with shaking. Since water takes too much time to reach room temperature, place a tray (for example of polystyrene) on the lid of the water bath to isolate the Plastic Trays from the hot water and use the agitation function of the water bath only.

- Add **1 ml** of **Conjugate**.
- Incubate **15 minutes** while shaking (about 50 rpm) at room temperature.
- Remove the solution, add **1 ml** of **Wash solution B** and rinse briefly.
- Remove the solution and add **1 ml** of **Wash solution B**.
- Incubate for **5 min** while shaking (about 50 rpm) at room temperature.
- Remove the solution and add **1 ml** of **Wash solution B**.
- Incubate for **5 min** while shaking (about 50 rpm) at room temperature.
- Remove the solution and add **1 ml** of **Color developer**.
- Incubate for **15 min** in the dark, while shaking (about 50 rpm) at room temperature.
- Remove color developer and wash strips with distilled water.
- Let strips dry in the dark on absorbent paper before proceeding with the interpretation of results.

⇒ AUTOMATIC DETECTION

Detection can be performed by using automatic instrument (*Tecan ProfiBlot 48 with software ver. 1.xx*): select “**Tipizza P**” program. In this case **10 µl PCR product + 10 µl DNAT** and **1,5 ml** of each reagent should be used for each sample.

⇒ RESULTS

Validation


A positive reaction of the uppermost Conjugate Control (CC) line indicates the correct function of the detection step. This line has to be always present and its intensity should be the same among different strips.

The colour intensity among the other positive lines on the strips may vary. This doesn't affect results interpretation.

Interpretation of results

- Align the CC line of the strip with the CC line on the transparent interpretation card. This line should always stain positive.
- The Amplification Control (AC) lines (lines I and II) contain universal probes for conserved sequences of the 5'UTR and Core region respectively. Their positive reaction ensures the addition of the amplified material and indicates the amplification of the corresponding region.
- Read the number of all positive lines for each strip on the Transparent interpretation card: the corresponding HCV typing result can be determined by using the Interpretation Chart (read also the section on diagnostic sensitivity for additional information on genotypes).

Strip scheme:

The part A of the interpretation table shows the possible result combinations on the basis of the 5'UTR region; the part B shows the possible result combinations obtainable on the basis of the Core region.

- If the 5'UTR region of the sample falls in the combinations of genotypes **2,3,4,5, 6a-b or 7a** don't continue with the analysis of the Core region. It is possible to observe positivity on lines **14-18** of the Core, but these **must not** be taken into account in the results interpretation.
- If the 5'UTR region of the sample falls in the combinations of genotypes **1 or 6c-q**, the analysis has to be continued with the interpretation of the results on the bands 14-18 (Core) and report only the results associated with the type indicated by the Core region (e.g. if a sample is genotype 1b in the 5'UTR region, but **1a in the Core region, it should be classified as 1a**).

In the presence of genotypes 2,3,4,5 and 7, signals on the probes of the Core region can be absent or present; in this case they should not be interpreted in any way:

Core probes	Genotypes that can show signals
14	4, 5a
15	/
16	2c, 4a, 5a
17	3
18	2, (3), 5a

TROUBLESHOOTING

1. False negative or too weak signals

- Weak signals on the strip could be the consequence of:
 - Too high DNA concentration during the PCR (try to dilute it) especially for samples with a very high titer ($>8 \times 10^6$)
 - Too low DNA concentration during the PCR for samples with a very low titer ($< 5 \times 10^3$)
- If all the strips show weak signals it could be due to too high hybridization and/or Stringent washing steps water-bath temperature. **Check carefully the water-bath temperature.**
- Hybridization or Stringent Wash Solution were not pre-warmed or mixed carefully. In this case the remaining solution should be discarded and a new one should be used.

2. False positive signals

- Unspecific signals can occur if temperature during hybridization and/or Stringent washing steps is too low. **Check carefully the water-bath temperature.**
- Stringent Wash Solution was not pre-warmed and mixed carefully. In this case the remaining buffer should be wasted and a new one should be used.
- If the same non-specific pattern, corresponding to a given genotype, is visible on all strips, a contamination occurred. Repeat the assay with new solutions. If contamination events occur it is important to discover their origin: extraction, RT and PCR reactions or contamination derived from other amplicons. The use of HCV positive and negative controls can be useful to find out the cause of contamination.

3. Single line non homogeneous staining

Shaking rate during detection steps is too low. Check that strips are submerged totally.

4. Pattern different from regular pattern on the typing list

- Typing probes of 2 different types are present: designate the sample as co-infection.
- Other Patterns: contact your distributor. In rare cases, combinations of bands are generated as not interpretable. This may be due to the heterogeneity of the HCV genome and to recombinant strains.

⇒ PERFORMANCE CHARACTERISTICS

Diagnostic specificity

The specificity of the Gen-C 2.0 assay was determined analysing 103 HCV negative plasma samples. Samples tested were negative for HCV antibody and RNA. The test was performed according to instructions for use. No false positive results were observed, therefore Gen-C 2.0 Test Specificity is 100%.

Diagnostic sensitivity

The diagnostic sensitivity of the kit Gen-C 2.0 was evaluated using plasma samples from HCV-positive subjects and amplified with kits cod NLM. AA910/48 ver. 2.0 . As reported in the interpretation chart, on the basis of the sequence analysis, the device is able to distinguish genotypes 1, 2, 3, 4, 6 and subtypes 1a, 1b, 2a/2c, 2b, 3a, 3b, 3c, 3k, 4a, 4b, 4c/4d, 4e, 4f, 4h, 5a, 6a/6b, 6g, 6f-q, 6m and 7a.

The device validation has been carried out with real samples of genotype: 1a, 1b, 2a/2c, 3, 3a, 4, 4c/4d, 5a and 6a. Viral loads were in a range of 2×10^3 and 1×10^7 IU/ml.

The Gen-C 2.0 assay ability to accurately genotype HCV positive samples was evaluated by comparing the results with those based on the previous version of the kit (for genotypes 2-5) and sequencing (for genotype 1 and 6) of individual samples; sequencing results were determined comparing them to HCV database and phylogenetic analysis. For genotypes 5a and 6a reference panels were used.

Analysis has been carried out at level of genotypes and at level of subtype for 1a and 1b. Of the 215 samples analyzed, all gave positive and correct results. The overall sensitivity for the assay is therefore 100%.

115/115 genotype 1 samples, were correctly subtyped with AC004 2.0 and were concordant with the sequencing method. Therefore the accuracy for Gen-C 2.0 assay for subtype 1a and 1b is 100%.

Summary of the analyzed samples and concordance:

Genotype	Nr of samples tested	Concordant to reference method (Gen-C, sequencing or panel)
1a	25	25
1b	90	90
2a/2c	37	37
3a	22	22
3	10	10
4c/4d	21	21
4	2	2
5a	6	6
6a/6b	2	2

Reproducibility

- **INTRA-ASSAY REPRODUCIBILITY (within series precision):**
Intra-assay reproducibility was evaluated using 3 samples of 3 different genotype in 3 replicates. The same 3 samples were tested with two different batches by two different operators.
- **INTER-ASSAY REPRODUCIBILITY (series to series precision):**
Inter-assay reproducibility was evaluated using 3 samples of different genotype by three different operators in different runs and days, using three different batches.

The test was performed according to the instructions for use, the precision was evaluated for all the steps of the procedure. The reproducibility for Gen-C 2.0 assay is 100%.

LOT

Lot
number

REF

Catalogue
number

Temperature
limitation

Use by

Consult instruction
for use

European
Conformity

Manufacturer

Irritant

IVD

In vitro
diagnostic device

Sufficient for

Keep away from
sunlight